

Specification Sheet

34998 Sigma-Aldrich

Acetonitrile

HPLC Plus, ≥99.9%

Synonym: ACN, Cyanomethane, Ethyl nitrile, Methyl cyanide

- CAS Number <u>75-05-8</u>
- Linear Formula CH₃CN
- Molecular Weight 41.05
- Beilstein/REAXYS Number 741857
- EC Number 200-835-2
- MDL number MFCD00001878
- eCl@ss 39031501
- PubChem Substance ID 329755104
- NACRES NA.03

CH₃CN

SKU-Pack Size	Availability	Pack Size	Price (SGD)	Quantity
34998-1L	Estimated to ship on 28.04.2021 - FROM	1 L	199.24	0
34998-2.5L	Estimated to ship on 05.05.2021 - FROM	2.5 L	297.56	0
34998-6X1L	Estimated to ship on 04.05.2021 - FROM	6 x 1 L	905.62	0
34998-4X2.5L	Estimated to ship on 04.05.2021 - FROM	4 x 2.5 L	991.87	0

Properties

Related Categories <u>Acetonitrile and Acetonitrile Solutions</u>, <u>Analytical</u>

Reagents, Analytical/Chromatography, Chromatography Reagents, HPLC Plus

Grade Solvents,



Sigma-Aldrich_®

Quality Level	<u>100</u>
grade	HPLC Plus
vapor density	1.41 (vs air)
vapor pressure	72.8 mmHg (20 °C)
assay	≥99.9%
form	liquid
autoignition temp.	973 °F
expl. lim.	16 %
baseline drift measuring range	≤12 mAU (210 nm)
application(s)	HPLC: suitable
	IR spectroscopy: suitable
impurities	≤0.0001% free alkali (as NH₃)
	≤0.0002% non-volatile matter
	≤0.001% free acid (as CH₃COOH)
	≤0.01% water (Karl Fischer)
	≤0.5 ppb fluorescence (quinine) at 254 nm
	≤0.5 ppb fluorescence (quinine) at 365 nm
halogenated residue	10 ng/L (as heptachlor epoxide)
refractive index	n20/D 1.344 (lit.)



Sigma-Aldrich.

bp	81-82 °C (lit.)	
mp	−45 °C (lit.)	
density	0.786 g/mL at 25 °C (lit.)	
HPLC-gradient	≤0.2 mAU at 254 nm	
	≤1 mAU at 210 nm	
λ	H₂O reference	
UV absorption	λ: 195 nm A _{max} : ≤0.10	
	λ: 200 nm A _{max} : ≤0.02	
	λ: 228 nm A _{max} : ≤0.005	
	λ: 235 nm A _{max} : ≤0.0044	
	λ: 250 nm A _{max} : ≤0.0044	
	λ: 400 nm A _{max} : ≤0.0044	
Featured Industry	Food and Beverages	
storage temp.	room temp	
SMILES string	CC#N	
InChl	1S/C2H3N/c1-2-3/h1H3	
InChI key	WEVYAHXRMPXWCK-UHFFFAOYSA-N	
	Show Fower Proportion	

Show Fewer Properties



Sigma-Aldrich_®

Description

General description

Acetonitrile (MeCN) is widely employed as a solvent in various studies and has high dielectric constant (37.5). Various purification procedures to obtain different grades of acetonitrile for use in different studies (polarography, spectroscopy, etc.) have been reported by many researchers.^[2]

Application

Acetonitrile has been used to compose the gradient system for the HPLC detection of diphenylmethane diisocyanate, a polymeric degradation product.^[1]

New Quality UPLC Plus Acetonitrile is suitable for all applications including:

Mass spectrometry, HPLC, TLC, GPC, UV, IR, GC, Pesticide Residue Analysis, and RNA/DNA Synthesis and Purification.

Packaging

1, 6×1, 2.5, 4×2.5, 4, 4×4 L in glass bottle

20 L in Nowpak™

20, 50 L in Pure-Pac™ 2

200 L in Pure-Pac™ 1

View returnable container options.

200 L in Pure-Pac™ 2

Preparation Note

Product filtered through a 0.2 µm filter

Other Notes

Pure-Pac® II containers require the Micromatic MacroValve coupler for dispensing solvents, **Z560723**.

Important notice

- The article number 34998-4X2.5L will be discontinued. Please order the single bottle 34998-2.5L which is physically identical with the same exact specifications.
- The article number 34998-6X1L will be discontinued. Please order the single bottle 34998-1L which is physically identical with the same exact specifications.

Recommended products

Discover LiChropur reagents ideal for HPLC or LC-MS analysis

Legal Information

Pure-Pac is a registered trademark of Sigma-Aldrich Co. LLC

